

EUROSPIRO SP-01, SP-11

Bluetooth & USB Spirometer

Spirometro Bluetooth e USB



LUMED®
www.lumed.com



***SP-01 & SP-11 are compatible with BOFAP !
The only cardboard AVB filter for spirometry .
Designed, patented, and produced by LUMED.***

EUROSPIRO SP-11 & SP-01 is the new solution for spirometry according to BTPS. It features A Bluetooth or USB interface, bi-directional turbine, and a colour display.

EUROSPIRO SP-11 & SP-01 support forced, relaxed and MVV vital capacity measurements. 3 measurement modes. Real-time monitoring of device connection and correction of body temperature and saturated pressure (BTPS).

EUROSPIRO SP-11 & SP-01 support flow-volume curve (volume velocity loop), and volume-time curve (vital capacity-time curve). Over 50 analysis parameters such as FVC, PEF, FEV1, FEV1/FVC, FEV6, lung age, etc. Detection, automatic diagnosis and end-of-exhalation.

Main features of EUROSPIRO SP-01 & SP-11:

- * No licences required
- * BOFAP compatible
- * Color display (SP-11)
- * Powered by just two AA batteries, alkaline or rechargeable (SP-11)
- * Bluetooth (SP-11) or USB (SP-01) interface
- * Android app (SP-11)
- * Uses our 3 liters calibration syringe

***SP-01 e SP-11 sono compatibili con BOFAP !
L'unico filtro AVB in cartone per spirometria.
Ideato, brevettato e prodotto da LUMED.***

EUROSPIRO SP-11 e SP-01 è la miglior soluzione per la spirometria secondo BTPS. Dispone di interfaccia Bluetooth o USB , turbina bidirezionale e display a colori.

EUROSPIRO SP-11 e SP-01 supporta misure di capacità vitale forzata, rilassata e MVV. 3 modalità di misurazione. Supporta il monitoraggio in tempo reale della connessione del dispositivo e la correzione della temperatura corporea e della pressione satura (BTPS).

EUROSPIRO SP-11 e SP-01 supporta la curva flusso-volume (anello velocità-volume) e la curva volume-tempo (curva capacità vitale tempo). Rileva fino a 50 parametri di analisi come FVC, PEF, FEV1, FEV1/FVC, FEV6 età polmonare, ecc. Supporta il rilevamento, la diagnosi automatica e la fine dell'espiazione.

Principali funzioni di EUROSPIRO SP-11, SP-01:

- * Non serve nessun tipo di licenza
- * Compatibile con BOFAP
- * Display a colori (SP-11)
- * Alimentazione con due sole batterie AA, alcaline o ricaricabili (SP-11)
- * Connessione al PC Bluetooth (SP-11) o USB (SP-01)
- * Applicazione Android (SP-11)
- * Usa la nostra siringa di calibrazione



Via Serio 10 20073 Opera ITALIA
tel +39 02 57606751
info@lumed.com

Via Senio 36/40 47121 Forlì ITALIA
tel +39 0543 702380
ordini@lumed.com