

TSTR0nn/1nn Disposable Nose Clips "EUROSPIRO"



N.B. Product is subject to change without prior notice. The actual product might differ from the pictures here shown.

Description

Nose clips are recommended for spirometry tests in order to reduce the risk of introducing a technical error. Our rubber- and foam-padded nose clips help prevent patients from inhaling or exhaling through the nose during spirometry.

In contrast to manual occlusion, the use of a nose clip can enhance patient comfort and allow the patient to focus on FVC, MVV, VC, IC and other spirometry maneuvers.

The nose clips are made of white thermoplastic, latex-free and CE marked. The tension in the nose clip body keeps it firmly in place while the rubber and foam pads are designed for maximum patient comfort. The nose clips are disposable and intended for single-patient use only.

Blue Large Nose Clip

Blue nose clips are larger and wider than white nose clips and are available with foam or rubber pads and in quantities of 10 or 100 per case.

White Nose Clip

White nose clips are narrower and smaller than blue nose clips and are available with foam or rubber pads in quantities of 10 or 100 per case.

Technical features

LUMED reference	Nose clip pads material	Nose clip body material	Color	Dimensions	Sale packaging
TSTR001/10	Polyvinyl	Polyacetate	White	70 x 40 mm	10 pcs in cardboard box
TSTR001/100	Polyvinyl	Polyacetate	White	70 x 40 mm	100 pcs in cardboard box
TSTR006/10	Ethylene-vinyl acetate	Polyacetate	White	70 x 40 mm	10 pcs in cardboard box
TSTR006/100	Ethylene-vinyl acetate	Polyacetate	White	70 x 40 mm	100 pcs in cardboard box
TSTR070/10	Polyvinyl	Polyacetate	Blue	70 x 45 mm	10 pcs in cardboard box
TSTR070/100	Polyvinyl	Polyacetate	Blue	70 x 45 mm	100 pcs in cardboard box
TSTR090/10	Ethylene-vinyl acetate	Polyacetate	Blue	70 x 45 mm	10 pcs in cardboard box
TSTR090/100	Ethylene-vinyl acetate	Polyacetate	Blue	70 x 45 mm	100 pcs in cardboard box