

## EP-LU4003 Universal Cart "EUROCARD"



N.B. Product is subject to change without prior notice. The actual product might differ from the pictures here shown. Arm and wooden handle are optional accessories.

### Description

LUMED "EUROCARD" Universal Medical Cart is a product line for equipment. The 1.5 mm thick steel structure that guarantees durability and reliability. Its painted surfaces are resistant and very easy to clean. Several extra accessories are available. The arm for patient cable has two hinges, and it contains the patient cables completely, for a better protection and usability. A built-in accessory bay is present in the rear side. It is equipped with four 70 mm diameter wheels and, two of them, with brake. The EUROCARD is tested for stability on an inclined plane up to 12°. It is supplied as an assembly kit includes the necessary tools and instructions.

## Technical features

<b>Manufacturer</b>	LUMED srl Via Serio, 10 - 20073, Opera (MI) - Italy Via Senio, 36 - 47121, Forlì (FC) - Italy
<b>Wheels</b>	4 wheels 70 mm diameter (2 wheels with brakes)
<b>Shelves</b>	Trolley, simple version with base and 4 wheels (2 with brakes), shelf and rear panel with compartment
<b>Containers</b>	Built-in bay, in the posterior side
<b>Colour</b>	Single color light-grey RAL 7035
<b>Dimensions</b>	40 x 45 cm (L x W) 90 cm (total height from pavement)
<b>Weight</b>	15 kg
<b>Payload</b>	15 kg each shelf; 40 kg max capacity
<b>Extra accessories</b>	<ul style="list-style-type: none"> <li>- EP-LU4002 – Adjustable arm for patient cable (RAL 7035)</li> <li>- EP-LU4004 - Additional shelf (RAL 7035)</li> <li>- EP-LU4005 - Storage compartment with flap door (RAL 7035)</li> <li>- EP-LU4010 - Upper shelf frame (RAL 7035)</li> <li>- EP-LU4015 - Antistatic wheels</li> <li>- EP-LU4016 - Wooden handle (RAL 7035)</li> <li>- EP-LU4017 - Probe support frame (RAL 7035)</li> <li>- EP-LU4019 - White steel basket</li> </ul>